LIROMA

RED LIGHT THERAPY

Manual

Portable PRIME

- Red light, Infrared light (NIR)
- Six wavelengths

Language index

English: Page 3-17


Nederlands: Pagina 18-32

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Français: Page 48-61

RED LIGHT THERAPY

Table of contents

Introduction.....	4
Red light therapy.....	5
Specifications.....	6
Benefits of Liroma red light therapy.....	7
Instructions for use.....	8
Product structure.....	9
Control panel	10
Which setting should I choose	11
Pulse function.....	11
Brightness	12
Installing the handle	12-13
Warnings.....	14-16
Care and cleaning	16
About Liroma	17

Thank you for your purchase!

Thank you for choosing Portable PRIME red light therapy!

Liroma wants its customers to feel healthy, fit and relaxed. To achieve this, Liroma has a wide range of high quality red light therapy products. Our products help with complaints such as: fatigue, poor circulation, pain in muscles and joints, pain due to underlying diseases, poor sleep, skin problems, stress and many other ailments. Our products target both mental and physical ailments.

Note: Red light therapy lamps do not get hot! During treatment, however, you may feel heat and sweat inside.

Note: Near-infrared light is invisible to the naked eye. It may appear as if some LEDs are not working, but this is completely normal.

Your wellness journey starts here!

Settle in, enjoy the moment, and let the soothing glow of red light support your body and mind. Here's to feeling stronger, sleeping better, and living each day with more vitality.


With love,

LIROMA

Red light therapy

Red light therapy uses a combination of red light and near-infrared light (NIR) to stimulate and repair the body at a cellular level. This therapy works by allowing the light to penetrate deep into the skin, with red light targeting the upper skin layers and NIR light reaching deeper into muscles and joints. A key effect of this light therapy is the stimulation of mitochondria—the energy factories of our cells—helping them produce energy more efficiently in the form of ATP (adenosine triphosphate).

With this increase in cellular energy, cells can recover more quickly, leading to reduced inflammation and improved blood circulation. Red light also promotes collagen production, which can contribute to firmer, more youthful skin and improved skin texture. Thanks to the combination of red light and NIR light, users often experience benefits such as reduced muscle and joint pain, enhanced skin health, faster muscle recovery, and increased energy levels. These effects make red light therapy a popular choice for those looking to support their well-being naturally.

**SCAN and learn how red light
therapy works**

Specifications

Portable
PRIME

Model	Portable Prime
Number of LED's	12 dual-chip LEDs
Wavelength (nanometer)	610 / 630 / 660 / 810 / 850 / 1060
Light spectrum	Red light / Near-infrared light (NIR)
Size	150 × 95 × 23 mm
Weight	2 kg
Power consumption	8 W
Brightness levels	10% / 45% / 70% / 100%
Pulse frequency	10 Hz – 9999 Hz
Irradiance	141 mW/cm ² (at 7.5 cm)
Timer	10 minutes (default)
Battery capacity	10,000 mAh
Charging time	Approx. 5 hours
Operating time	Approx. 45 min per full charge
Material	Aluminum
Color	Silver

Benefits of Liroma red light therapy

Product benefits

- ✓ Combination of red and near-infrared (NIR) light
- ✓ 6 wavelengths: 610 / 630 / 660 / 810 / 850 / 1060 nm
- ✓ Compact and portable design
- ✓ User-friendly control panel
- ✓ Adjustable brightness levels
- ✓ Pulse mode
- ✓ Timer
- ✓ Rechargeable battery
- ✓ LED lamps (energy-efficient)


Health benefits

- ✓ Reduction of muscle and joint discomfort
- ✓ Support of muscle recovery after exercise
- ✓ Improvement of skin condition
- ✓ Anti-inflammatory effect
- ✓ Stimulation of blood circulation
- ✓ Support of wound healing
- ✓ Reduction of acne and other skin issues
- ✓ Increase in cellular energy production
- ✓ Relief from local tension and stiffness


Instructions for use

- **Ideal usage frequency:** Start with 3–4 sessions per week to allow your body to adapt. After 2–4 weeks, you may increase usage if desired.
- **Maximum distance:** Use the device at approximately 7–10 cm from the treatment area.
- **Recommended treatment example:** Standard treatment time is 10 minutes per session.
- Ensure that the body part being treated is fully exposed to the light from the lamp and is not covered by clothing.
- Always consult a doctor before use in case of contraindications, medical conditions, or if you are unsure.
- **Note:** The lamp does not emit heat! It is designed solely for LED therapy.

Product structure

Control panel

Power button

Turn the device on or pause the treatment.

Timer

The default timer is set to 10 minutes. Press and hold (3s) the power button to change the time.

Brightness

Adjust the brightness of the lights with this button. There are L1(10%), L2(45%), L3(70%), L4(100%) brightness level, to match different sensitive skin.

Battery indicator

Shows remaining battery power. This device has 10,000 mAh rechargeable battery. To fully charge the battery requires approximately 5 hours.

The device supports approximately 1 hour of operating time per full charge.

Pulse

Frequency adjustable from 10-9999Hz, which are proved to have positive performance on cells.

Mode

Choose between the following light modes: RED, NIR, RED + NIR

We recommend using both red and near-infrared light simultaneously for most purposes.

Which setting should I choose?

Red light (610, 630, 660nm)

- Improves skin health
- Stimulates blood circulation
- Stimulates collagen and elastin production

NIR (near-infrared light) (810, 850, 1060nm)

- Reduces inflammation
- Relieves muscle and joint pain
- Restores deep tissues such as muscles, joints, and organs
- Stimulates cellular energy in deeper tissues

We always recommend turning on all the lights at the same time for the best results.

PULSE-function

Pulsed light in red light therapy penetrates deeper into the tissue, stimulates cell repair, and increases ATP production. It enhances circulation, reduces heat stress, and can activate specific biological responses at different frequencies.

Additionally, research suggests it may offer neurological benefits, such as in Alzheimer's and Parkinson's diseases. A combination of continuous and pulsed light may yield the best results.

We recommend experimenting with pulsed settings between 1 Hz and 10 Hz to find what feels most comfortable for you.

Brightness

A dimmer function on a red light therapy panel offers users the flexibility to adjust the light intensity. This provides greater comfort, especially for sensitive skin types, and allows for customized therapies.

Additionally, a lower intensity helps during longer sessions and can reduce energy consumption, which in turn extends the lifespan of the device.

Installing the handle


Portable PRIME can be used handheld or with the included handle for more stable positioning.

Follow the steps below to install the handle correctly.

Step 1


Place the back of the light panel facing upward.

Position the handle so that it aligns with the designated mounting holes on the back of the device.

Step 2


Insert the screws through the handle and into the mounting holes.

Step 3

Tighten the screws securely until the handle is firmly attached.

Ensure the handle is properly fixed before using the device.

Warnings

Warnings and precautions

- Do not use this device on babies.
- Do not allow sleeping, disabled, or unconscious individuals to use this device.
- Do not use the device for more than 20 minutes per therapy session.
- Do not use the device on open wounds.
- Do not use the device in combination with liniments, ointments, or balms that contain heat-producing ingredients, as this may cause burns.
- Do not attempt to modify or repair the device yourself. If it is damaged or not functioning properly, please contact us as soon as possible at info@liroma.nl.

This device is electrical. To prevent electric shock and other hazards, please follow these instructions:

- Do not wash the product with water or other liquids.
- Do not place the device in water or in areas where it may be exposed to water or other liquids.
- Do not expose the device to running water.
- If the device comes into contact with or is submerged in water or any other liquid, immediately unplug it from the socket. If using the device near water, we recommend using a socket with a ground fault circuit interrupter (GFCI).
- Do not touch the device or its plugs with wet hands.
- Do not use or touch the device with wet or damp skin.
- Do not lift the device by the power cord.
- Do not use the device if the power cord is damaged.
- Do not use the device on or near heated surfaces.
- Do not use the device in areas where it may be exposed to flammable or combustible products or fumes, as this could cause a fire or explosion.
- Do not store the device at temperatures above 140°C (284°F).

Additional warnings

- Do not use the device if it is damaged, as continued use of a damaged device may cause electric shock or injury.
- Do not use extension cords unless they are designed to handle the total power of the device.
- Do not unplug the device by pulling on the cord. Always grasp the plug itself when removing it from the outlet.
- This device is not suitable for individuals with physical, neuro-sensory, or mental impairments, or those lacking experience and/or knowledge.
- Do not use the lamp in combination with implanted electronic devices, such as artificial hearts or artificial lungs, without first consulting your doctor.
- Do not use the lamp on pregnant women, individuals with severe heart conditions, critically ill patients, or emergency cases without consulting your doctor.
- Discontinue treatment immediately if you experience hypersensitivity or allergic skin reactions due to medication, cosmetics, or food.
- Under no circumstances should you open the device or attempt to repair it yourself, as this may cause injury and void your warranty. If the device is damaged or not functioning properly, contact the retailer or info@liroma.nl as soon as possible.

Do not use this product without first consulting your doctor if you suffer from:

- Severe illnesses such as heart disease, acute inflammation, thrombosis, blood clotting disorders, adrenal suppression, systemic lupus erythematosus, or malignant conditions.
- Hypersensitivity to infrared light or have edema.
- A condition related to nerve demyelination, such as carpal tunnel syndrome, multiple sclerosis, or ulnar neuropathy.
- Pregnancy.
- Certain medications or plant extracts that cause photosensitivity are not compatible with light therapy. (Examples: antibiotics, antidepressants, Roaccutane, MTX, antifungal medication, and St. John's Wort).

Important notes

- Red light is visible, while NIR infrared light is invisible to the naked eye. However, with NIR infrared light, you may perceive a slight glow or dimming, which is completely normal.

Care and cleaning

Follow the instructions below to clean and maintain the device properly and safely.

- Switch off the device before cleaning
- Unplug the device from the charging cable
- Remove the device from any accessories before cleaning
- Use a dry or slightly damp microfiber cloth
- Do not use chemical cleaning agents or alcohol
- Do not spray liquid directly onto the device
- Be careful when cleaning around LED openings and vents
- Do not expose the device to water

LIROMA

RED LIGHT THERAPY

Liroma has been a specialist in red light therapy for years. We offer lamps that provide full-body exposure, as well as smaller devices that allow for targeted treatment on specific body parts. Additionally, we have red light devices specifically designed for the back and face. Visit our website to explore all our products.

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