
Instruction Manual

LIROMA MANUAL STAND

For Horizontal and Vertical Use

Language Index

English: Page 01-15

Nederlands: Pagina 16-28

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Français: Page 42-54

Table of Contents

Stand components & hardware	06
Product overview	07
Assembly instructions	08
Using the RTL stand	12
Care and maintenance	14
Warranty	15

Welcome to Liroma Red Light Therapy

Thank You for Choosing Liroma Red Light Therapy!

At Liroma, we are committed to helping our customers feel healthier, fitter, and more relaxed. To support this goal, we offer a wide range of high-quality red light therapy products designed to promote overall wellbeing.

The RTL Stand is designed to provide a stable, adjustable mounting solution for your compatible Liroma red light therapy panel, allowing you to position it comfortably over a bed, sofa, or other suitable location.

The stand features adjustable height and angle settings, making it easy to position your therapy panel for a wide range of treatment applications while improving comfort and ease of use.

The RTL Stand is compatible with all Liroma RTL therapy panels. It is not compatible with the Liroma Ultra, Multi or products from other manufacturers.

For safe operation, use this stand only with compatible Liroma red light therapy panels and do not exceed the maximum supported load specified in this manual.

Important Notes:

- The RTL Stand is designed to improve the positioning and mobility of your compatible Liroma therapy panel, allowing it to be adjusted to different heights and angles for comfortable use.
- When repositioning the stand, lower the therapy panel to a comfortable height to improve stability and maneuverability.
- When the stand is not in use, store it in a clean, dry indoor location to help protect its finish and moving components.
- Periodically inspect the stand for signs of wear and clean it with a soft, dry or slightly damp cloth. Avoid using abrasive cleaners or solvents.

Thank you for trusting Liroma as part of your wellness journey.

With appreciation,

The Liroma Team


WARNINGS:

Improper assembly or installation may result in serious personal injury or property damage. Read and follow all instructions before assembling or using this product.


- If you do not understand these instructions, or if any parts are damaged or missing, contact Liroma Customer Support before assembly.
- Do not assemble or use this product if any component is damaged or missing. Use only genuine replacement parts supplied by Liroma.
- The maximum supported load is 9.8 kg (21.6 lbs.). Do not exceed this weight limit.
- This stand contains moving and adjustable parts. Keep hands and fingers clear of moving components during adjustment to avoid injury.
- Use this stand only for its intended purpose and only with compatible Liroma red light therapy panels.
- Do not modify this product or use non-approved accessories. Liroma is not responsible for damage or injury resulting from product modifications, improper assembly, incorrect installation, or exceeding the maximum weight capacity.
- **WARNING:** Children must not climb, hang, or play on the stand or mounted therapy panel. The stand may tip over, resulting in serious injury or death.
- Do not place toys, remote controls, or other attractive objects on the mounted therapy panel, as this may encourage children to climb on the stand.
- Use only with compatible Liroma therapy panels. Installing incompatible products or products exceeding the maximum weight capacity may cause the stand to become unstable, resulting in injury or property damage.

Note: The illustrations in this manual are provided for reference only. Your product may vary slightly in appearance. Follow the assembly instructions to ensure correct installation.

Stand Components & Hardware

A Long Screw
× 4

B Medium
Screw × 4


C Short Screw
× 5


D Hand
Tightening
Knob × 9

E Left Base Tube
× 1


F Right Base
Tube × 1

G Support Rod
× 1


H Lift Rod
× 1

I Base Connecting
Tube × 1

J Rotating Rod
× 1

K Mounting
Hook × 2

L Swivel
Caster
× 4


M Protective Cap
× 2

Instruction
Manual × 1

Product Overview


1. Panel Mount Lock Knob
2. Mounting Rail
3. Mounting Hook
4. Tilt Lock Knob
5. Base Assembly
6. Pivot Lock Knob
7. Height adjustment knob
8. Support Rod
9. Swivel Caster

Assembly Instructions


Step 1: Assemble the Base

1. Align the Left Base Tube (E) and Right Base Tube (F) with the Base Connecting Tube (I).
2. Insert the Medium Screws (B) through the mounting holes and tighten all the screws securely using the appropriate tool.


Step 2: Install the Casters

1. Turn the assembled base upside down.
2. Screw the Swivel Casters (L) into the threaded holes at the ends of the base tubes. Tighten each caster securely.

Step 3: Install the Support Rod

1. Insert the Support Rod (G) into the centre opening of the base.
2. Secure the support rod using Long Screws (A).
3. Turn the base over and insert Short Screw (C) into the bottom of the base to lock the support rod in place. Tighten all screws securely.

Step 4: Install the Lift Rod

1. Insert the Lift Rod (H) into the Support Rod (G). Adjust the lift rod to the desired height by aligning the adjustment holes.
2. Insert the Protective Caps (M) into the top openings of the support rod and lift rod.
3. Secure the lift rod using Hand Tightening Knobs (D). Tighten both knobs firmly to lock the lift rod at the selected height.


Note:

- Align the adjustment holes before installing the hand-tightening knobs.
- Ensure both knobs are fully tightened before mounting the therapy panel.

Step 5: Install the Rotating Rod

1. Insert the Rotating Rod (J) into the top of the Lift Rod (H).
2. Secure the rotating rod using Hand Tightening Knob (D) on the side of the lift rod.
3. Tighten the knob firmly to lock the rotating rod in position.


Mounting Orientation

The rotating rod can be installed in two orientations depending on how you intend to mount the therapy panel.


- **Horizontal Panel Mounting**

Position the rotating rod so that the row of mounting holes faces upward. This orientation allows the panel to be mounted horizontally.

- **Vertical Panel Mounting**


Rotate the rod 90° so that the row of mounting holes faces sideways. This orientation allows the panel to be mounted vertically.

Important: Verify that the mounting holes face the correct direction before tightening the hand-tightening knob. This will ensure the therapy panel can be installed correctly in the desired orientation.

Step 6: Mount the Therapy Panel

1. Attach the Mounting Hooks (K) to the back of the therapy panel using Short Screws (C). Tighten the screws securely.
2. Lift the therapy panel and carefully place the mounting hooks onto the Rotating Rod (J).
3. Align the mounting holes on the hooks with the holes in the rotating rod.
4. Secure the panel using Hand Tightening Knobs (D). Tighten both knobs firmly.
5. Check that the therapy panel is level and securely fastened before use.


Important: Before each use, ensure both hand tightening knobs are fully tightened and the therapy panel is securely mounted. Do not use the stand if the panel is loose or unstable.

Using the RTL Stand

The RTL Stand allows your compatible Liroma therapy panel to be positioned in different orientations to suit your treatment needs.


Vertical Orientation:

1. Use the panel in the vertical position when treating larger areas of the body, such as:
 - Back
 - Spine
 - Legs
 - Full-body standing treatments
2. Adjust the height and tilt of the stand so the treatment area is centered in front of the therapy panel.

Horizontal Orientation:

1. Rotate the mounting assembly so the panel is in the horizontal position when treating areas while lying down or when a wider treatment area is preferred, such as:
 - Chest
 - Abdomen
 - Shoulders
 - Full-body treatments over a bed or sofa
2. Adjust the height and angle to position the panel comfortably above the treatment area.

Care and Maintenance

Proper care will help maintain the stability and appearance of your Liroma Stand.

- Wipe the stand with a soft, dry or slightly damp cloth. Do not use abrasive cleaners, solvents, or harsh chemicals.
- Keep the stand dry at all times. Do not immerse any part of the stand in water or allow liquids to enter the adjustment mechanisms.
- Regularly inspect all screws, hand-tightening knobs, and mounting components. Tighten any loose hardware before use.
- Check the casters periodically to ensure they rotate freely and remain securely attached. Remove any dirt or debris that may affect movement.
- Inspect the stand for signs of wear, damage, or deformation. Do not use the stand if any component is damaged or unstable.
- When moving the stand, push it slowly on a level surface. Avoid sudden impacts, uneven surfaces, or obstacles that could cause the stand to tip.
- Store the stand indoors in a clean, dry location away from excessive moisture, direct sunlight, and corrosive environments.
- Use only genuine Liroma replacement parts. Do not modify the stand or replace components with non-approved parts.

Important: Always ensure the therapy panel is securely mounted and all hand-tightening knobs are fully tightened before each use.

Warranty

A limited two-year warranty from the date of purchase covers the Liroma Stand. This warranty applies only to the original purchaser and is non-transferable.

This warranty covers defects in materials and workmanship under normal use during the applicable warranty period. If a defect arises during the warranty period, Liroma will, at its discretion, repair or replace the defective product or defective part at no additional charge.

If the product cannot be repaired or replaced, Liroma may provide an equivalent replacement product or refund the original purchase price, subject to applicable consumer laws.

For the first 180 days after the date of purchase, Liroma will cover the shipping costs associated with an approved warranty claim. After the first 180 days, the customer is responsible for the cost of shipping the product to Liroma for warranty service. Any repaired or replacement product will be covered for the remainder of the original warranty period.

This warranty does not cover damage resulting from:

- Misuse, abuse, or negligence.
- Improper installation or maintenance.
- Unauthorized modification or repair.
- Improper storage or transportation.
- Normal wear and tear.
- Theft, loss, or accidental damage.

Any attempt to disassemble, modify, or repair the product by anyone other than an authorized service technician may result in injury and will void this warranty. Liroma is not liable for any damage, injury, or loss resulting from unauthorized disassembly, repair, or attempted repair of the product.

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